



... Please, fill in completely and keep a copy ...

www.medimatch.ie

MediMatch dental laboratory IRELAND T: 01 830 9375

Lyon House Business Centre, Unit 6, 50 Nore Road, Dublin Industrial Estate, Glasnevin Dublin 11, D11 R978 • questions@medimatch.ie

CROWN & BRIDGE PRIVATE

for lab use only, real app.

..... - @
..... - @

For lab use only, Opened by:

1 | 2

1 Dr. _____ (name) Acc./PO no.: _____
Post code: _____ - _____ Mobile: _____

2 *No patient personal data*
Your case reference: _____ Age: _____ Initials: ____ / ____

3 **Delivery date** Always give us **10 full** working days
(d/m/y) _____ / _____ / _____ = 1 working day before the real appointment

4 **Margin (for porcelain bonded)**
 No metal margin
 Metal margin lingual / palatal (= standard)
 Metal all around by _____ mm (standard 0.2mm anterior, 0.5mm posterior)
 Metal backing
 Metal backing + metal palatal cusp
 Porcelain facing only

For lab use only, enclosed / missing:

I have disinfected the impression



! with: _____
by: _____

Pictures? E-mail to: questions@medimatch.ie

5 * See also dedicated order form

Restoration
 Full-ceramic*
 E.max
 Press (= very aesthetic) €140
Faster delivery or CAD-CAM products available at a surcharge
 Zirconia
 full contour (posterior only) €100
 Multi Layer (with stain & glaze) €115
 with porcelain (=aesthetic) €140
 Porcelain bonded
 non precious (Co-Cr)
 semi precious (Pd)
 precious, Gold (Au)
 Full Metal
 Co-Cr Pd Au
 Composite
 Post and core
 Co-Cr Pd Au
 Integral Separate


6 **Shade** - use VITA guide
 Bridge B
 Crown C
 Inlay I
 Onlay O
 Veneer V
 Wings on nr: _____



7 **Instructions** Porcelain **butt margin** on: _____
 Telescopic crowns on: _____
 Special **attachments** on: _____
 Wax up on: _____

PRIVATE only

Zirconia (Cad-Cam) is available as the original cercon zirconia - € 135



Return special delivery min. € 5
 This case is a remake case
 I have enclosed new or old components _____

Custom Made Device | Supplied in an unsterilized state | Terms and conditions apply | Visit www.medimatch.ie

Dental appliance information and delivery note: Dental product designed to satisfy the information, properties and detail of what has been prescribed by the above dentist. Purely for use for the patient described with the above reference. The product meets the requirements entailed in the Annex I and Annex VIII of the Medical Devices Directive and the Republic of Ireland Medical Devices Regulations act SI 252/1994. This product may have been produced in either or all three Ireland, UK or China by MediMatch. **Instructions for use, storing & handling:** It is highly recommended that the product is stored in a clean and safe environment if not used immediately. It is also advised that there must be no contact with materials, liquids or acids that could cause disfiguration or damage to the product. The product should not be subjected to extreme heat. Where applicable, you should take care not to damage the dental piece(s) when removing from its model. THIS DENTAL APPLIANCE IS SUPPLIED IN AN UNSTERILIZED STATE